AT51500

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Thank you for purchasing our product. Please read this instruction manual carefully before use. Do not operate the product before you have read and understood this manual. If you do not understand any of the warnings, cautions, or instructions, contact your healthcare professional or the distributor to avoid damage to the product or personal injury.

CAUTION: Check all components for damage caused during transportation. If any damage is found, DO NOT use the product. Contact the distributor for further information.

INTENDED USER GROUP: This product is intended for individuals affected by diseases, dysfunctions, or injuries that require treatment, rehabilitation, or compensation through the use of this device (see the section “Intended Purpose” in this manual). The product may be purchased directly by the user or based on the recommendation of a physician, therapist, or other healthcare professional. Regardless of how the product is obtained, factors such as available sizes, required functions, variants, indications and contraindications for use, as well as the manufacturer’s information, should always be taken into account.

INTENDED PURPOSE: The product is designed to assist visual perception in individuals with reduced visual acuity. It enables magnification of the observed object, making reading, identifying details, performing precision tasks, and examining small components easier. The product can be used in home, outpatient, or hospital environments as an optical aid for everyday use or for diagnostic activities that do not require contact with the patient’s body.

INDICATIONS FOR USE: The magnifying glass is particularly recommended for persons:
with reduced visual acuity caused by refractive errors (e.g., presbyopia, myopia);
with ocular diseases leading to limited vision (e.g., cataract, glaucoma, age-related macular degeneration – AMD);
undergoing vision rehabilitation or recovering after ophthalmic surgery;
requiring visual support for precision activities such as reading fine print, viewing medication labels, inspecting skin, surgical sutures, or operating medical/prosthetic devices; for educational or diagnostic purposes – as an auxiliary tool for medical personnel or caregivers.

CONTRAINDICATIONS

Use of the magnifying glass is not recommended in the following cases:
physical or mental limitations (e.g., severe visual impairment) that prevent safe handling of the product;
hypersensitivity to bright light (photophobia), especially when using models with LED lighting;
impaired spatial orientation or visual-motor coordination that could hinder safe use;
cases where an ophthalmologist advises against optical aids (e.g., after certain eye surgeries or in the treatment of photosensitive glaucoma);
environments with direct exposure to sunlight, due to the risk of concentrated light causing burns or fire.

ADDITIONAL NOTES: The magnifying glass is not a substitute for corrective lenses or ophthalmic treatment. If visual difficulties persist, consult an ophthalmologist. The product is intended for external observation only – do not use it for intraoral, ophthalmic, or dermatoscopic examinations.

LIST OF COMPONENTS: Magnifying Glass, Instruction manual

PRODUCT FEATURES

The lens is made of PMMA and utilizes Fresnel lens technology, providing a natural, clear image and high light transmittance.

The ergonomic handle is coated with soft rubber for a comfortable, non-slip grip. The light source consists of 32 cool-white LEDs and 16 warm-white LEDs, evenly distributed in four rows around the lens. They provide bright, even illumination with the option to switch between light modes (cool, warm, or mixed). A large battery compartment holds 3 AAA (1.5 V) batteries for extended use. The rotating handle allows flexible positioning and easy folding for transport and storage. The hinge mechanism includes a spring-loaded locking system made of wear-resistant material to ensure durability.

TECHNICAL SPECIFICATIONS

Materials: ABS frame, rubber handle, PMMA lens, Fresnel optical system

Lens dimensions: 18.5 × 11 cm

Folded dimensions: 20.8 × 16.5 × 2.5 cm

Unfolded dimensions: 30 × 13.5 × 2.5 cm

Magnification: 5×

Power supply: 3 × AAA (1.5 V)

Weight: 133 g

Light source: 32 cool LEDs + 16 warm LEDs

OPERATION


Hold the frame with one hand and adjust the handle with the other to the desired angle. Place the object on a stable surface, bring the magnifier closer, and slowly move it along the visual line until the image becomes sharp. If needed, press the LED button to turn on the light. You may select lighting based on ambient conditions to protect your eyesight.

LIGHT MODES

Press once – Cool white light/ Press twice – Warm white light/ Press three times – Mixed light (cool + warm)/ Press four times – Turn off

BATTERY REPLACEMENT

Open the battery compartment cover. Remove used batteries and insert 3 new AAA batteries, ensuring correct polarity. Close the cover securely.

CAUTION – SERIOUS INCIDENT REPORTING

If a serious incident related to the device has occurred, which has led, might have led, or could lead to any of the following:

- a) death of a patient, user, or other person;
- b) temporary or permanent serious deterioration in the health of a patient, user, or other person;
- c) a serious public health threat;

the incident must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident.

ADDITIONAL WARNING: If pain, allergic reactions, or any unusual or unclear symptoms occur during use, discontinue use and consult a healthcare professional. It is forbidden to use the product in any way other than its intended purpose.

CLEANING: Clean the lens surface with a dry microfiber cloth. If necessary, wipe gently with a damp cloth using a mild detergent, then clean water. Do not use alcohol, gasoline, or other aggressive chemicals.

STORAGE: Store the product in a dry place. Do not expose it to direct sunlight or high temperatures.

DISPOSAL: After the end of its service life, dispose of the product as regular household waste. Observe local waste segregation regulations. Protect the environment.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o., 03-068 Warsaw, ul. Zawisłńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.
10. The warranty does not cover parts and components subject to natural wear and tear, such as wheels, tires, rubber pads, fuses, etc., the wear of which results from normal use of the product. The warranty period for batteries is 6 months.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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