AT53056, AT53057, AT53058

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PURPOSE/ INDICATIONS

Product is used to strengthen the lumbar spine. Sewn-in reinforcements as well as additional side supports provide additional stiffening. The elastic, „breathable“ material allows the natural warmth of the body to be maintained and improves blood circulation. The spinal orthosis is recommended for a variety of spinal disorders, e.g. lumbar discopathy; degenerative spine disease; root syndromes (sciatica).

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling.

WARNING

This product is intended for personal use only. Do not use the product if there is any damage on the supplied accessories, components or on the product itself. Do not wear the orthosis if ointments or warming fluids have been applied to the area. If a rash or irritation occurs, remove the orthosis and discontinue use. Consult a doctor immediately. Use the product as directed by your doctor.

INSTRUCTIONS FOR USE

The product may be available in different sizes. Please refer to the size chart to find the correct size. Take the product out of the package, then put on the lumbar spine and fasten the Velcro fastener. The possibility of stretching the elastic bands with a velcro end will ensure proper pressure adjustment. Unless otherwise advised by your doctor, wear the device for no more than 4 hours a day.

WASHING INSTRUCTIONS

Hand wash with mild soap at up to 40 °C. Dry laid flat, away from heat and direct sunlight. Do not bleach, iron or dry-clean. Do not tumble dry.

AT53057, AT53058 - remove stiffeners before washing.

STORAGE

Store in a dry and cool place, away from direct sunlight.

SPECIFICATIONS

AT53056	Waist circumference (cm)	Height (cm)	Length (cm)
M	70-93	23	98
L	90-103	23	108
XL	100-113	23	118
XXL	110-123	23	128
XXXL	120-133	23	138

AT53057	Waist circumference (cm)	Height (cm)	Length (cm)
M	70 – 85 cm	26	90
L	80 – 95 cm	26	100
XL	90 – 105 cm	26	120
XXL	100 – 115 cm	26	120
XXXL	100 – 125 cm	26	130

AT53058	Waist circumference (cm)	Height (cm)	Length (cm)
M	60-80	32	88
L	77-93	32	98
XL	87-110	32	108
XXL	107-123	32	118
XXXL	117-130	32	128

NOTES:

In the event of a device-related „serious incident“ that has directly or indirectly led to, could have led to, or is likely to lead to any of the following:

- (a) death of the patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

the above „serious incident“ must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTE: In the event of pain, allergic reactions or other distressing symptoms, unclear to the user, associated with the use of the medical device, consult a healthcare professional.

NOTE: It is forbidden to use the product other than in accordance with its intended use

ATTENTION: The manufacturer is not responsible for any damage caused by negligence in cleaning, storage or as a result of not following the recommendations contained in these instructions for use.

HOW TO DISPOSE OF THE PRODUCT AT END-OF-LIFE

Once the device has been taken out of service, the medical device can be disposed of as normal municipal waste.

SCOPE OF SUPPLY: Orthosis, instructions for use

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in this instruction manual). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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